
X-Ray Tube Housing Assembly

Gaine radiogène
Röntgenröhrengehäuse
Coraza del tubo de rayos X
X射线管组件

453575069095 Rev A

DA 200 ULTRA Series
Série DA 200 ULTRA
DA 200 ULTRA Reihe
Serie DA 200 ULTRA
DA 200 ULTRA系列

DUNLEE

INSTRUCTIONS FOR USE

Section 1.0:

Safety Information, Informations relatives à la sécurité, Sicherheitsinformationen, Información de seguridad, 安全资讯..... p. 4-23

Section 1.1:

Intended Use, L'emploi destine, Bestimmungsgemäße Verwendung, Intención de uso, 使用用途 p. 6

Section 1.2

General, Note relative à la sécurité, Sicherheitshinweis, Advertencias de seguridad, 常规..... p. 7-11

Section 1.3

Radiation Protection, Radioprotection, Strahlenschutz, Protección contra radiaciones, 辐射防护..... p. 12

Section 1.4:

Electrical Safety, Sécurité électrique, Elektrische Sicherheit, Seguridad del sistema eléctrico, 电气安全 p. 13

Section 1.5

Thermal Safety, Sécurité thermique, Thermische Sicherheit, Seguridad térmica, 热安全 p. 14

Section 1.6:

Electromagnetic Compatibility (EMC), Compatibilité électromagnétique, Elektromagnetische Verträglichkeit (EMV), Compatibilidad electromagnética (CEM), 电磁兼容性(EMC)..... p. 15-16

Section 1.7:

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Copyright

© 2012 DUNLEE

All rights reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright-holder. DUNLEE reserves the right to make changes in specifications or to discontinue any product, at any time without notice or obligation, and is not liable for any consequences resulting from the use of this publication.

Equipment is subject to change without notice. All changes will be in compliance with regulations governing manufacture of medical equipment. Document originally drafted in English.

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1.1 Intended Use

This X-ray Tube Housing Assembly is designed and manufactured to be compatible with specific diagnostic imaging equipment. Your local X-ray Tube supplier can provide confirmation on request,

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1.2 Safety Notice

This manual is designed to make it possible for you to work with the x-ray tube housing assembly in a safe manner. You may operate the X-ray tube housing assembly only in compliance with the safety instructions in this manual and not use it for purposes other than for which it is intended. The X-ray system may only be operated by persons who have the necessary expertise in radiation protection and who have been instructed in how to operate the X-ray system.

It is always the user who is responsible for compliance with the regulations applying to installation and operation of the X-ray unit

- You must never use the x-ray tube housing assembly or any associated diagnostic imaging equipment if it has any electrical, mechanical or radiological defects. This particularly applies to faulty indicators, displays, warnings and alarms.
- If the assembler wishes to connect the X-ray tube housing assembly to other equipment, components or assemblies and if it is not apparent from the technical data whether it can be safely combined with such equipment, components or assemblies, the user must ensure that the safety of the patient, operating staff, third parties and the environment is not affected by the planned combination by consulting the manufacturers involved or by making inquiries from an expert.
- Dunlee and Philips are responsible for the safety features of its products only if persons explicitly authorized to do so by Dunlee or Philips have performed maintenance, repairs and modifications.
- As with any technical appliance, this equipment requires
 - correct operation
 - regular, competent maintenance
 - care
- If you operate the X-ray equipment, and hence the x-ray tube housing assembly, incorrectly or if the user fails to maintain it properly, Dunlee and/or Philips cannot be held responsible for any malfunctions, damages or injuries.
- Safe operation of the X-ray tube assembly is only guaranteed when it is used according to its specification. In case the specification limits are disregarded, there is the danger of oil leaks and expelled component parts as a result of implosion/explosion. In such cases the manufacturer of the X-ray tube assembly excludes any liability. Any guarantee claims for this product are rejected.
- The safety circuit of the X-ray tube assembly that prevents the switching ON of radiation when the specified temperature limit of the X-ray tube assembly is exceeded shall be connected before the X-ray tube assembly is set to work for the first time. It shall neither be removed nor modified.

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1.3 Radiation Protection

- Make certain before every X-ray exposure that all necessary radiation precautions have been taken. You can find information about radiation precautions in the instructions for use for the individual X-ray systems with which you are using this X-ray tube assembly.

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1.4 Electrical Safety

- Only trained maintenance staff may remove the covers from the X-ray tube assembly.
- This X-ray tube assembly may only be used in medical rooms which meet the requirements of relevant national and international standards and laws.
- You must not use this X-ray tube assembly in areas where there is danger of explosion.
- Detergents and disinfectants, including those used on the patient, can create explosive gas mixtures.
- Please observe the relevant regulations.

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1.5 Thermal Safety

To prevent damage due to thermal overloading, precautions must be taken to ensure that the X-ray tube assembly is not operated outside its specified load parameters. This is the only way to eliminate all risk to the patient, operating staff, third parties and the environment.

The X-ray tube assembly, consisting of the X-ray tube and X-ray tube housing, is one component in the X-ray system. The power supply and the frame rate are generated and controlled by the X-ray system. The user information of some X-ray generators indicate the thermal state of the X-ray tube.

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1.6 Electromagnetic compatibility (EMC)

In accordance with its intended use, this electronic component is tested as part of a diagnostic imaging system according to the laws governing EMC, which defines the permitted emission levels from electronic equipment and its required immunity against electromagnetic fields.

Electronic apparatus that satisfies the EMC requirements is designed so that under normal conditions there is no risk of malfunction caused by electromagnetic interference. However, in the case of radio signals from high-frequency transmitters with a relatively high transmitting power, the risk of electromagnetic incompatibility when operated in close proximity to electronic apparatus cannot be totally ruled out.

In unusual circumstances unintended functions of the apparatus could be initiated, possibly giving rise to undesirable risks for the patient or user. For this reason, all kinds of transmission with mobile radio equipment should be avoided. This also applies when the apparatus is in STANDBY mode. Mobile telephones must be switched OFF in designated problem zones.

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1.7 Electromagnetic emission

The responsibility to guarantee compliance with the electromagnetic emission limit values and the fulfillment of all associated regulations and laws lies with the company responsible for the installation of the X-ray tube assembly.

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1.8 Disposal

 Take-back, proper disposal and recovery of the Medical Device takes place in accordance with the European WEEE Directive (Waste Electrical and Electronic Equipment) and the requirements of national legislation.

Our state-of-the art x-ray equipment is manufactured according to safety and environmental protection procedures. Under normal operation, there are no risks to persons or the environment.

To comply with regulations it is necessary to use materials which may be harmful to the environment and therefore have to be disposed of in a proper manner.

This X-ray tube assembly may contain materials which are toxic. For this reason, you must not dispose of the X-ray equipment together with industrial or domestic waste.

The tube housing assembly is filled with insulation oil. The disposal of this oil must be done in accordance with your local regulations. Refer to the disposal section in this manual for details. SHELL Diala AX does NOT contain Dioxin.

 This x-ray tube housing assembly may contain beryllium as part of the x-ray tube window material. Disposal of products with beryllium must be through appropriate channels. The Product Specification Sheet in section 4 will indicate if beryllium is present.

Dunlee/Philips supports you in the disposal of the X-ray tube assembly in accordance with valid regulations and laws.

We strive to be an environmentally conscious with the products and the material returned for disposal. Consequently, Dunlee carefully disassembles and processes parts and materials according to earth friendly objectives. Some of these materials may be used in the manufacture of new parts and assemblies as part of this Green initiative . All Dunlee products will achieve the same high level of quality and functionality our customers have come to expect.

Contact your Dunlee/Philips organization for proper handling and disposal.

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1.10 Transportation

Extreme care must be taken when transporting x-ray tube housing assemblies. It is recommended to transport x-ray tube housing assemblies in the original packaging material. Transportation using other packaging material or in any other way, may result in product damage and may void the warranty.

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1.11 Temperature and Pressure Safety Switch

The x-ray tube housing assembly is equipped with a thermal housing safety switch. Activation of the switch will inhibit additional exposures.

When this occurs, the unit will remain inoperable only until a safe operating temperature is restored. Depending on the cooling system, it may require 5 to 10 minutes cooling time to reduce the temperature.

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2.1 Installation

Only trained and qualified service personnel must perform installation of the x-ray tube housing assembly. The installation must be performed in accordance with the system documentation provided by the Original Equipment Manufacturer's installation, test and operating procedures. Consult the appropriate documentation for the removal and installation procedure.

This x-ray tube housing assembly is intended to be installed on a diagnostic imaging system that is fitted with an appropriate audible warnings and/or visual displays, where practicable, while radiation is being omitted.

Failure to follow the OEM procedure for removal and installation of the x-ray tube housing assembly will void the warranty.

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2.2 X-ray Tube Housing Conditioning

It is imperative for proper operation and extended life that the x-ray tube housing assembly be properly conditioned after a period of non-use. This procedure is defined in the Original Equipment Manufacturer's operator's manual. Failure to follow this procedure may reduce the life of the x-ray tube housing assembly and may void the warranty terms.

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2.3 General Maintenance

The responsibility to recommend preventive maintenance in accordance with all applicable regulations and laws lies with the company responsible for the installation of the X-ray tube assembly.

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2.4 Corrective Maintenance

Corrective maintenance on X-ray tube assemblies is allowed to be performed only by the manufacturer of the X-ray tube assembly.

Faulty components of the X-ray equipment which affect the safety of the X-ray tube assembly must be replaced by genuine spare parts.

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2.5 Planned Maintenance – General

X-ray tube assemblies contain mechanical components, which are subject to normal wear due to operation.

The correct setting of the electromechanical and electronic assemblies affects the functioning, image quality, electrical safety and exposure of the patient and medical personnel to radiation.

As with any technical appliance these X-ray tube assemblies also require:

- Regular checks by the operator.
- Regular planned and corrective maintenance.

By taking these precautions you maintain the operability and operational reliability of the system. As the user of the X-ray unit, you are obliged according to accident prevention regulations, the medical products law and other regulations to perform such precautionary actions.

Maintenance consists of tests that the user can perform and maintenance that is performed under service agreements or by persons explicitly authorized to do so.

 X-ray tube assemblies contain mechanical components which are subjected to wear and tear due to operation. The correct setting of the electromechanical and electronic assemblies safeguards the functioning, image quality, electrical safety and exposure of the patients, operating staff, third parties and the environment.

The user must check the x-ray equipment for apparent defects (see section 2.7 checks performed by the operator). If operational defects or other departures from normal operational behavior occur, the x-ray unit must be switched off and the appropriate service organization informed. Operation of the x-ray equipment may only be resumed when repairs are completed. Operation using faulty components may lead to an increased safety risk or unnecessarily high exposure to radiation.

It is recommended that you perform the tests indicated in the table on a regular basis and have the unit serviced by authorized representatives at least once a year. For heavily used equipment, preventative maintenance should be scheduled more often.

These precautionary measures prevent personal injury and guarantee that the operator fulfills all commitments.

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2.6 Maintenance by Operator

 Apply the following measures of maintenance only if the X-ray tube assembly is directly accessible. The configuration of the X-ray equipment must not be changed.

Cleaning

Detergents with a high alcohol content make the material dull or cause it to crack.

 ***Switch OFF the mains voltage of the X-ray equipment before cleaning. Never use any caustic, solvent or abrasive detergents or polishes.***

Heed during cleaning:

- Be careful that no water or other fluids can get inside the X-ray equipment. This prevents short circuits in the electrical installations and corrosion on the components.
- Clean enamelled parts and aluminium surfaces only with a moist cloth and a mild detergent and wipe with a dry woollen cloth.

Disinfection

The disinfection method used must conform to valid legal requirements and guidelines for disinfection and protection against explosion.

Switch OFF the mains voltage before disinfection of the X-ray equipment. If you use disinfectants which create explosive gas mixtures. these must first have dispelled before restarting the X-ray equipment.

Wipe disinfection

Disinfect all parts of the X-ray equipment, including accessories and connecting cables, by wiping only with a soft cloth.

Spray disinfection

Spray disinfection is not recommended because disinfectant could get inside the X-ray equipment.

Disinfection by atomizing

Cover the X-ray equipment carefully with a sheeting when it has cooled down. After the disinfection mist has subsided, you may remove the sheeting and disinfect the X-ray equipment by wiping with a soft cloth.

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3.1 Compatibility

The x-ray tube housing assembly described in this manual meets the provisions of the Medical Directive 93/42 EEC (93) and FDA 21 CFR 1020.30.

If you have further questions regarding the applicable national or international standards, please address them to:

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Patents

The x-ray tube housing assembly is manufactured under one or more of the following United States patents:

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3.3 Carrier Damage

In the unlikely event of damage due to handling by the carrier, it is important to follow a few specific steps in order to receive proper credit:

1. Inspect the x-ray tube housing assembly immediately upon receipt. Check both the packing and the product for physical damage.
2. If there is physical damage, call the carrier immediately and order a "Joint Inspection" of both the packing and the product.

Within the United States:

3. If the x-ray tube housing assembly has been shipped freight pre-paid, return the product in accordance with the instruction on the Return Form along with a copy of the Joint Inspection Report to Dunlee.
4. If the x-ray tube housing assembly has been shipped via the customer's carrier, the customer must initiate the claim process directly with the carrier.

Outside the United States:

5. Please consult your local sales and service office for instruction for return.

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3.4 Operational Warning

If there is evidence of oil leakage from the x-ray tube housing assembly, heat exchanger, or hoses, immediately cease operation of the system and inform your service organization or supplier.

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3.5 Warranty

Warranties vary according to the specific tube type and manufacturer. Contact your local representative for detailed information. Written copies are available on request.

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3.6 Return Process

Use the re-usable packaging of the new delivered X-ray tube housing assembly for the return shipment of the defective X-ray tube housing assembly. This guarantees a transport of the defective X-ray tube housing assembly in a packing which has been developed for this purpose.

Complete the Tube Return Form:

If a paper form is supplied with the new X-ray tube assembly, follow the instruction on the Return Form, fill in all required information and send it together with the defective X-ray tube housing assembly to the designated facility. You may elect to provide the defective tube information via electronic format by completing the form found on the Dunlee.com webpage.

If only an electronic form is supplied on a CD or a memory stick, follow the instructions provided and send the memory stick together with the defective X-ray tube housing assembly to the designated facility.

Outside the United States:

If returning a product from outside the United States, please contact your local sales or service representative for return instructions.

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Philips Healthcare (Suzhou) Co.,Ltd